

Haemostatic Powder HemoCeeTM PLUS



- 100% plant based
- Hydrophilic particles
- 100% absorbable
- Off-the-shelf deployment
- Diverse range of applications

BioCer

HaemoCer™ PLUS- Derived from nature.

4 Step Application Technique “**DAPI**”

Dry- Remove excess blood from the wound site

Apply- Immediately apply HaemoCer™ PLUS

Press- Apply compression for 2 minutes

Irrigate- Irrigate well prior to removal of gauze



Reference Number	Product	Content
HFP202	HaemoCer™ PLUS HAEMOSTATIC POWDER 2g incl. Applicator	5 pcs
HFP205	HaemoCer™ PLUS HAEMOSTATIC POWDER 5g incl. Applicator	5 pcs
HFZ101	HaemoCer™ UNIVERSAL APPLICATOR 440 mm	5 pcs

BioCer

BioCer Entwicklungs-GmbH

Ludwig-Thoma-Straße 36c · 95447 Bayreuth · Germany

Tel. +49 (0921) 78 77 70-0 · Fax +49 (0921) 78 77 70 79

customerservice@biocer-gmbh.de · www.biocer-gmbh.de



Distributed by

 **surgical
specialties**

Cutting Edge Innovation

3/17 Rodborough Rd, Frenchs Forest NSW 2086
PO Box 6052, Frenchs Forest DC NSW 1640, Australia
Customer Service T: 1300 665 884 F: 1300 665 886
www.surgicalspecialties.com.au