

non-narcotic pain relief means

Better outcomes. Satisfied patients.

Product Ordering Information

Every day without



costs you more.



A KIMBERLY-CLARK HEALTH CARE COMPANY



As your trusted provider of market-leading, elastomeric ON-Q® pumps for more than 10 years, I-Flow® has expanded its portfolio to provide you with a comprehensive suite of products to meet your needs. ON-Q® pumps provide significantly better pain relief than narcotics alone.

ON-Q® Pain Relief System

ON-Q® PUMP KITS WITH ON-Q® SOAKER® CATHETER

Delivering local anesthetic automatically and continuously to the surgical site for days following surgery—providing effective pain relief.

Delivery Time† (in days)	Fill Volume x Flow Rate	Soaker® Catheter 1 in (2.5 cm)	Soaker® Catheter 2.5 in (6.5 cm)	Soaker® Catheter 5 in (12.5 cm)	Per Package
2	100 ml x 2 ml/hr	PS2502	PS6502	PS12503	1
5	270 ml x 2 ml/hr	—	PS6503	—	1
2¼	270 ml x 5 ml/hr	PS2504	PS6504	PS12504	1
3	270 ml x 4 ml/hr (dual; 2 ml/site)	—	—	—	1
4	100 ml x 1 ml/hr	—	—	—	1
3½	400 ml x 5 ml/hr	—	PS6506	PS12506	1
1½	400 ml x 10 ml/hr	—	PS6507	PS12507	1
4	400 ml x 4 ml/hr (dual; 2 ml/site)	—	—	PS12508	1
2	100 ml x 2 ml/hr (dual; 1 ml/site)	PM039	—	—	1

ON-Q® PUMP KITS WITH ON-Q® SILVERSOAKER® CATHETER

The same ON-Q® pumps you know and trust with the added benefit of an antimicrobial silver coated catheter.

Delivery Time† (in days)	Fill Volume x Flow Rate	Soaker® Catheter 2.5 in (6.5 cm)	Soaker® Catheter 5 in (12.5 cm)	Soaker® Catheter 10 in (25 cm)	Per Package
2	100 ml x 2 ml/hr	PS6502A	—	—	1
5	270 ml x 2 ml/hr	PS6503A	—	—	1
2¼	270 ml x 5 ml/hr	PS6504A	PS12504A		1
3	270 ml x 4 ml/hr (dual; 2 ml/site)	PS6505A	PS12505A	PS25005A	1
4	100 ml x 1 ml/hr	—	—	—	1
3½	400 ml x 5 ml/hr	—	PS12506A	PS25006A	1
4	400 ml x 4 ml/hr (dual; 2 ml/site)	—	PS12508A	PS25008A	1

ON-Q® PUMPS

ONDEMAND® bolus button is an added feature that allows patients to give a 5 ml bolus device in addition to the continuous flow rate.

Delivery Time† (in days)	Fill Volume x Flow Rate	Fenestrated Catheter 1 in (2.5 cm)	Per Package
Varies	270 ml x 2 ml/hr (5 ml bolus) (60 minute refill)	PMB01	1

ON-Q® PUMPS

The Select-A-Flow® controller allows you to adjust the flow rate according to your patient's needs.

Delivery Time† (in days)	Fill Volume x Flow Rate	Fenestrated Catheter 1 in (2.5 cm)	Per Package
Varies	Select-A-Flow 270ml x 1-7ml/hr, 6.5cm	PSAF02	1

PUMPS FOR CONTINUOUS PERIPHERAL NERVE BLOCKS

ON-Q® pumps for CPNB slowly infuse local anesthetic near a nerve. ONDEMAND® bolus button is an added feature that allows patients to give a 5 ml bolus device in addition to the continuous flow rate.

The Select-A-Flow® controller allows you to adjust the flow rate according to your patient's needs.

Delivery Time† (in days)	Fill Volume x Flow Rate	Product Code	Per Package
2	400 ml x 10 ml/hr	CB002	1
ONDEMAND®			
< 3	400 ml x 5 ml/hr with 5 ml bolus (60 minute refill)	CB003	1
Select-A-Flow®			
Varies	400 ml x 2-14 ml/hr	CB004	1
Varies	600 ml (dual; 1-7 ml/site)	CB6007	1
ONDEMAND® & Select-A-Flow®			
Varies	400 ml x 2-14 ml/hr with 5 ml bolus (30 minute refill)	CB006	1

Surgical Products

ON-Q* SOAKER* & ON-Q* SILVERSOAKER* CATHETER EXPANSION KITS

Patented design provides an even distribution of local anesthetic for optimal post-op pain relief.

Fill Volume x Flow Rate	Soaker* Product Code	SilverSoaker* Product Code	Per Package
Soaker 2.5 in (6.5 cm), T-Peel Needle 3.25 in	PM010	PM010-A	1
Soaker 5 in (12.5 cm), T-Peel Needle 6 in	PM020	PM020-A	1
Soaker 1 in (2.5 cm), T-Peel Needle 3.25 in	PM030	PM030-A	1
Soaker 10 in (25 cm), T-Peel Needle 8 in	—	PM040-A	1
Soaker 7.5 in (19 cm), T-Peel Needle 8 in	—	PM050-A	1

ON-Q* TUNNELERS

ON-Q* tunnelers provide flexibility to enhance optimal catheter placement.

Description	Product Code	Per Package
Tunneler Sheath 11 GA x 8 in	ACC07S	10
Tunneler Sheath 17 GA x 8 in	ACC08S	10
Tunneler Sheath 11 GA x 12 in	ACC09S	10
Tunneler Sheath 16 GA x 12 in	ACC10S	10
Disposable Tunneler & Sheaths 17 GA x 8 in	T17X8	5
Disposable Tunneler & Sheaths 11 GA x 8 in	T11X8	5
Disposable Tunneler & Sheaths 11 GA x 12 in	T11X12	5
Disposable Fluid Tunneler & Sheaths 11 GA x 8 in	FT11X8	5

ON-Q* ACCESSORIES

Description	Product Code	Per Package
T-Peel Needle 6 in	ACC02	10

There are inherent risks in all medical devices. Please refer to the product labeling for **Indications**, **Cautions**, **Warnings**, and **Contraindications**. Failure to follow the product labeling could directly impact patient safety. Physician is responsible for prescribing and administering medications per instructions provided by the drug manufacturer. Refer to www.iflo.com for product safety Technical Bulletins.



KIMBERLY-CLARK HEALTH CARE COMPANY

Rx only

*Registered Trademark or Trademark of Kimberly-Clark Worldwide, Inc.

© 2013 KCWW. All rights reserved.

SSIFLOWPI 042014

Distributed by



Cutting Edge Innovation

4B Ride Way, PO Box 301-218

Albany, Auckland, New Zealand

Customer Service T: 0800 665 884 F: 09 447 1685

www.surgicalsurgicalspecialties.co.nz

3/17 Rodborough Rd, Frenchs Forest NSW 2086

PO Box 6052, Frenchs Forest DC NSW 1640, Australia

Customer Service T: 1300 665 884 F: 1300 665 886

www.surgicalsurgicalspecialties.com.au