

# PediNail™

Pediatric & Adolescent Femoral Nail

## Your patient. Your nail. Your solution.

### ***Anatomically designed***

to meet the unique needs of pediatric and adolescent patients

### ***Lateral trochanteric entry point***

to facilitate a safe and efficient procedure

### ***9mm proximal diameter***

allows for easy insertion and reduced reaming

### ***Nail diameters from 7-10mm***

for a broad patient offering



## Features and Benefits

- Complex 3 dimensional geometry, results in anatomically appropriate design
- Advanced surgical technique reduces the danger to blood vessels near the piriformis fossa
- Smaller diameter allows for easier insertion without the need for excessive reaming
- Minimally invasive technique facilitates early mobility and return to normal activities
- Complete set of easy to use instrumentation

The Adverse Effects and Contra-Indications listed in the IFU are important considerations. These Risks and general surgical Risks should be explained to the patient prior to surgery.

## PediNail™ Pediatric and Adolescent Intermedullary Nail\*

### 3 Proximal locking hole options

- One Oblique hole (greater-to-lesser trochanter)
- One Transverse (9 & 10mm nail only)
- One Recon hole (into femoral neck)

### Available in Multiple Diameters & Lengths – Right & Left

- Diameters: 7, 8, 9 & 10mm
- Lengths: 20 - 38 cm
- All nails are cannulated



Proximal lateral bend accommodates lateral trochanteric entry

Anterior bow to mimic normal pediatric anatomy for age 8-16 years

Distal lateral bend to ease insertion lateral to the greater trochanter, minimizing medial cortex interference

2 medial to lateral distal locking holes

### PediNail™ Screws

4.0/4.5 Fully-Threaded Cortical Screws



4.5 Recon Screws



(Available in multiple lengths)

1 anterior to posterior distal locking hole

\*Patent Pending



OrthoPediatrics, Children Are Not Just Small Adults, PediPlates, PediBlade, PediLoc, PediFlex, PLEO, PediNail, Scwire and the Pedi logo are trademarks of OrthoPediatrics Corp.

OrthoPediatrics, Children Are Not Just Small Adults, PediPlates, PediLoc, Scwire and the Pedi logo are registered trademarks in the United States.

2850 Frontier Drive • Warsaw, IN 46582 • ph: 574.268.6379 or 877.268.6339 • fax: 574.268.6302 • [www.OrthoPediatrics.com](http://www.OrthoPediatrics.com)

OrthoPediatrics Corp. ©2012 Part# SS-1500-01-01 Rev D